



On-site Training at  
Takeo Provincial Hospital  
Laboratory  
November 12-15, 2018

# Quality Control Overview

Ong Siew Kim  
MS, MAACB, FAACB, MBA, PBT(ASCPi), EdD  
Project Coordinator

# Quality Control

Establish range for new lot:

- ✓ At least 20 values
- ✓ At least 2 weeks or 10 working days, and
- ✓ Preferably over at least 4 weeks or 20 working days

# QC range calculation

1. Calculate **mean** = 20 data results new QC lot.
2. Use your lab **historical CV** from old QC.
3. Calculate SD.

$$\text{(Historical) CV} = \frac{\text{SD} * 100}{\text{Mean}}$$

$$\text{SD} = \text{Mean} * \text{CV}\%$$

4. Set “temporary” QC range (new lot).
5. Set up LJ chart.

# Practical 1

## QC Handling process

- Preparation
- Storage
- Thawing

# Handling process- Preparation

1. Take out QC bottle from a fridge.
2. Leave on table.
3. Take diluent or water place on the table.
4. Wait for 15 mins. Set timer.
5. Use pipette to re-constitute the QC.
6. Leave for 15 mins. Set timer.
7. Roll gently 15 mins. Set timer.
8. Ready for use.



# Handling process- Storage

1. Prepare your aliquot tubes.
2. Calculate the number of aliquots.
3. Keep in freezer (master log for reagent)
4. Details on storage box: Labelling
  - a. Content and quantity, concentration or titer (manufacture lot number, level)
  - b. Storage requirements (temperature)
  - c. Date prepared
  - d. Expiry date
  - e. Med Tech Initial
5. Fill in QC log book.
6. Fill in QC LJ chart.
7. Limit work to <1 hour.

## Handling process- Thawing

1. Take out a frozen aliquot from a freezer
2. Leave on table.
3. Wait for 15 mins. Set timer.
4. Mix gently. Use a micropipette.
5. Limited to room temperature exposure to <1 hour



Questions?